ACCU-CHEK, ACCU-CHEK RAPID-D LINK, ACCU-CHEK SPIRIT, ACCU-CHEK SPIRIT COMBO, COMBO and AKKY-CHEK are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2023 Roche Diabetes Care

Made in Poland

Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.accu-chek.com

Last update: 2023-08


04537220001/0

ACCU-CHEK® Rapid-D Link











Instructions for use

Infusion set


ACCU-CHEK®

Suitable infusion sites / Geeignete Infusionsstellen / Sites de perfusion appropriés / Siti d'infusione idonei / Zonas de infusión adecuadas / Locais de infusão adequados / Geschikte infusieplaatsen / Passende infusionssteder / Egnede infusionssteder / Lämpliga infusionsställen / Sopivat infuusiokohdat / Κατάλληλα σημεία έγχυσης / Odpowiednie miejsca infuzji / Vhodná místa infuze / Vhodné miesta infúzie / מקומות מתאימים להחדרת עירוי / Infúzió beadásra alkalmas helyek / Primerna mesta za infuzijo / Locuri de administrare a perfuziei adecvate / Подходящие места установки / Sobivad infusioonikohad / Piemērotas dūriena vietas / Tinkamos infuzijos vietos / مواضع التسريب المناسبة

Instructions for use	6	en
Gebrauchsanweisung	16	de
Instructions d'utilisation	26	fr
Istruzioni per l'uso	36	it
Instrucciones de uso	46	es
Instruções de utilização	56	pt
Gebruiksaanwijzing	66	nl
Brugsanvisning	76	da
Bruksanvisning	86	no
Bruksanvisning	94	sv
Käyttöohje	104	fi
Οδηγίες χρήσης	112	el
Instrukcja obsługi	122	pl

Návod k použití	132	cs
Návod na použitie	140	sk
הוראות שימוש	157	he
Használati útmutató	158	hu
Navodila za uporabo	168	sl
Instrucțiuni de utilizare	176	ro
Инструкция по использованию	186	ru
Kasutusjuhend	198	et
Lietošanas pamācība	206	lv
Naudojimo taisyklės	216	lt
إرشادات الاستخدام	237	ar

PRECAUTION

Risk of serious harm

The disregard of the instructions for use can prevent the product from working as intended. This can lead to various harms, including hypoglycemia, hyperglycemia, injury, or infection.

Carefully read and follow the instructions for use. Also, read the User's Manual of your insulin pump.

Intended Use

The infusion set is intended for the subcutaneous infusion of insulin delivered by an insulin pump.

Indication

For use in the therapy of diabetes mellitus.

Intended User

People with diabetes and healthcare professionals or caregivers who assist people with diabetes in using the product.

Target Group (applicable for)

Infants, children, adolescents, and adults with diabetes.

Contraindications

There is no clinical experience for the use

with neonates.

Clinical Benefits

The infusion set enables pump therapy when it is connected to a compatible insulin pump. The infusion set enables connected insulin pumps to fulfill their purpose, which is to provide an intensive insulin therapy.

Limitations

The device is intended for single use only. The device is intended to be used only with insulin or insulin analogs.

Additional Materials Required

The infusion set is intended to be used together with the following products only:

Insulin pumps

- the Accu-Chek Spirit insulin pump
- the Accu-Chek Spirit Combo insulin pump

Cartridges

- the Accu-Chek Spirit 3.15 ml cartridge system
- the Insuman® Infusat 100 IU/ml solution for injection in a cartridge (3.15 ml, Sanofi Aventis)

Device Description

Cannula

- A** Tubing
- B** Connector
- C** Steel cannula
- D** Needle cover
- E** Protective film
- F** Adhesive pad
- G** Handling aid

Transfer set

- H** Tubing
- I** Luer-lock connector
- J** Protective cap
- K** Connector

Protective cap

- L** Protective cap

The infusion set includes a cannula and a transfer set. The cannula is made of steel and it is pointed and sharp. The cannula is inserted at a 90° angle into the subcutaneous tissue. A short tubing is part of the cannula.

The transfer set must be connected to the cannula and through a standard luer-lock connector to the insulin pump.

The cannula, the transfer set, and the protective cap (L) are sterile, single use products that have been sterilized with ethylene oxide.

Contents of the Pack

There are 3 packages with different contents. The cannula ([REF](#) 04541006, 04541014, 04541022, 04541049), the transfer set ([REF](#) 04541057, 04541065, 04541073, 04541081), and the protective cap (L) ([REF](#) 04540999) are packed in their own separate packages.

General Safety Information

A  **WARNING** indicates a foreseeable serious hazard.

A  **PRECAUTION** describes a measure you should take to use the product safely and effectively or to prevent damage to the product.

WARNING

Risk of suffocation

This product contains small parts that can be swallowed.

Keep the small parts away from small children and people who might swallow small parts.

Risk of strangulation

If the tubing of the transfer set is wrapped around the neck, the tubing can strangulate you.

Pay attention to the tubing. The transfer set is available in different lengths. Use the shortest possible length.

Risk of hyperglycemia

If the connections in the infusion set are not tight, the insulin delivery can be insufficient.

If you feel unwell, if you suspect insufficient insulin delivery, or if you unexpectedly smell insulin, check whether your blood glucose value increases. Also, check the infusion set for leaks.

Risk of hyperglycemia or injury

If the cannula is used in areas with electromagnetic radiation, the cannula can heat up and become occluded or the cannula can shift in the subcutaneous tissue. This can lead to hyperglycemia or injury.

Do not use the cannula in areas with electromagnetic radiation, for example, during a MRI scan. Stop and disconnect your pump. Pull the cannula out of the

skin.

Risk of infection

- Expired products can be unsterile.
Use only products that are within the use by date. The use by date is printed next to the  symbol on the packaging.
- If the sterile packaging is damaged or already open, the product can be unsterile and lead to infection.
Visually inspect the sterile packaging. If the packaging is damaged or already open, do not use the product. The sterile packaging is marked with the  symbol.
- If a cannula was already used, it can transmit infections.
Use a cannula only once.

PRECAUTION

Risk of hyperglycemia, increased pain, or injury

If you insert the cannula into a body part that is not suitable as infusion site, the insulin delivery can be insufficient or increased pain, or injury can occur. Select a suitable infusion site. Ask your healthcare professional where you can insert the cannula.

Risk of hyperglycemia

- Several reasons can lead to an insufficient insulin delivery. After attaching or replacing the cannula or the transfer set, check your blood glucose level at least once within 1 to 3 hours. If your blood glucose level is unexpectedly high, follow the instructions in these instructions for use.
- Air in the infusion set can lead to an insufficient insulin delivery. Follow the instructions in these instructions for use on how to fill the infusion set completely or how to remove any air bubbles.

Risk of hyperglycemia or inflammation

Mechanical stress, for example, pressure or force, on the infusion set or the infusion site can lead to an insufficient insulin delivery or inflammation. Avoid mechanical stress on the infusion set and on the infusion site.

Risk of inflammation

In rare cases, the product can cause unexpected body reactions. If you experience skin irritation or redness, contact your healthcare professional.

Risk of infection, inflammation, or scarring

If you use infusion set components for a longer period than intended, infection, inflammation, or scarring can occur. Replace the cannula at least every 2 days. Replace the transfer set at least every 6 days.

If your blood glucose level unexpectedly increases or if an occlusion alarm occurs, check the infusion set for occlusion or leaks. If you are not sure, whether the infusion set is working properly, replace the infusion set immediately. If your infusion site is inflamed, replace the infusion set immediately. Select a different infusion site during the healing process.

If the head set or the transfer set was dropped or exposed to other mechanical stress, check it for damages, for example, cracks. If it is damaged, do not use it. Check the connection between the luer-lock connector (I) and the insulin pump adapter every 3 hours during the day and before you go to bed. Check your blood glucose level at least 4 times a day.

Check at regular intervals that the adhesive pad does not peel off from the infusion site and that the cannula is fully inserted. If the cannula is not fully inserted anymore, replace the cannula.

After you replace the infusion set or the cartridge and before you insert the cannula into the skin, check whether the infusion set is completely filled with insulin.

Suitable Infusion Sites

Select an infusion site with sufficient subcutaneous tissue, for example, on the abdomen or the back of your upper arm. See the illustration on page 2. The infusion site must be at least 5 cm away from the navel and the last infusion site used and must not be over a bone. The infusion site must not be exposed to rubbing or bumping, for example, through a waistband, belt, or safety belt. The infusion site must not have skin changes, for example, scars, injuries, bruises, rashes, or tattoos.

Rotate between different infusion sites each time you replace the cannula to avoid tissue changes from insulin. Think of

a pattern to select new infusion sites.

Attaching the Infusion Set

When you use the infusion set for the first time, your healthcare professional must be present.

- Wash your hands.
- Disinfect the chosen infusion site as instructed by your healthcare professional.
- Let your skin dry completely.

PRECAUTION

Risk of infection

If you do not disinfect the infusion site, an infection can occur.

Disinfect the selected infusion site as instructed by your healthcare professional before you insert the cannula into the skin.

See illustrations on page 1.

- 1** Hold both connectors (B and K) so that they form a cross. Plug the connectors into each other. Turn both connectors a quarter turn. You can only turn the connectors in one direction. For a tight connection, make sure that you hear and feel the

connectors click into place.
If you do not hear and feel both connectors click into place, replace the cannula and the transfer set.

- Unscrew the protective cap (J) from the luer-lock connector (I) of the transfer set.

2 PRECAUTION

Risk of hypoglycemia, hyperglycemia, or infection

If the luer-lock connector is damaged, the insulin delivery can be insufficient or uncontrolled or germs can enter the infusion set.

To avoid damaging the luer-lock connector, do not use any tools to screw the luer-lock connector into the insulin pump adapter. Only screw by hand.

Screw the luer-lock connector (I) tight into the insulin pump adapter.

3 WARNING

Risk of hypoglycemia

If you fill the infusion set while the infusion set is connected to your body, an uncontrolled insulin amount

can be delivered.

Never fill the infusion set while the infusion set is connected to your body.

Fill the infusion set as described in the User's Manual of the insulin pump until insulin drips from the steel cannula (C). You must fill up the infusion set (tubings and cannula) completely before you insert the cannula into the skin. The infusion set must not contain any air bubbles. Remove any air bubbles while filling the infusion set.

Fill volume of the infusion set (U100 insulin)

Length	Units	Volume
30 cm / 12"	approx. 7 U	approx. 70 µL
60 cm / 24"	approx. 10 U	approx. 100 µL
80 cm / 31"	approx. 13 U	approx. 130 µL
110 cm / 43"	approx. 18 U	approx. 180 µL

- 4 Hold the cannula by the handling aid (G).
Pull off both parts of the protective

film (E) from the adhesive pad (F).

- 5** Unscrew the needle cover (D) to remove it.
Do not touch the cannula to avoid injury and contamination.
- 6** Form a skin fold at the chosen infusion site.
Insert the steel cannula (C) quickly and at a 90° angle into the skin fold.
- 7** Press the adhesive pad (F) firmly onto the skin.
Check that the adhesive pad is properly attached to the skin.
- 8** Hold the adhesive pad (F) in place.
Bend the handling aid (G) to the side until it detaches itself.
- 9** Tape the excess tubing as a loop to the skin with an adhesive pad. This can prevent you from getting caught with the tubing and thus pulling on the steel cannula.

Disconnecting and Reconnecting the Insulin Pump

This infusion set allows you to disconnect your insulin pump, for example, for showering and swimming. You can reconnect the insulin pump without having

to replace the infusion set.

Consult your healthcare professional on how to compensate for any missed insulin while you are disconnected from the insulin pump.

Disconnecting the insulin pump

- Turn both connectors a quarter turn in opposite directions until they form a cross. You can only turn the connectors in one direction.
- Pull the connectors apart.
- Place a protective cap (L) on each connector.

Reconnecting the insulin pump

- Remove the protective caps from the connectors.
- Connect the tubing (H) and the cannula using the connectors as described in step **1**.

Replacing the Cannula

- Put your pump into STOP mode.
- Disconnect the pump.
- Carefully pull off the adhesive pad (F), working round the outside edge into the center.
- Pull the steel cannula (C) out of the skin at a 90° angle and completely pull off

the adhesive pad simultaneously.

- Connect the pump to a new cannula as described in step **1**.
- Then continue with step **3**.

Replacing the Transfer Set

- Put your pump into STOP mode.
- Disconnect the pump.
- Unscrew the luer-lock connector (I) of the transfer set from the insulin pump adapter.
- Connect a new transfer set to the pump as described in step **2**.
- Fill the new transfer set as described in the User's Manual of the insulin pump.
- Connect the new transfer set to the cannula as described in step **1**.

Discarding Used Parts of the Infusion Set

PRECAUTION

Risk of infection

A used product that has come into contact with human body fluids can transmit infections.

Discard all used parts of the infusion set as infectious material according to local regulations. For more information see the following instructions.

Discard the pointed steel cannula in such a way that it does not cause injury to anyone. For example, place the cannula into a sharps container.

For information on how to discard the product correctly, contact your local council or authority.

Operating and Storage Conditions

Use the infusion set between 5 and 40 °C. Store the infusion set in the unopened blister at room temperature. Temporarily, the temperature can be between 5 and 45 °C. Keep the product dry and away from sunlight.

Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

Customer Support

Australia

Roche Diabetes Care Australia Pty. Limited
Insulin Pump Support: 1800 633 457
australia.insulinpumps@roche.com
www.accu-chek.com.au

Roche Diabetes Care South Africa (Pty) Ltd.
Hertford Office Park, 90 Bekker Road
Vorna Valley, 1686

South Africa

Tel: +27 (11) 504 4600
Email: info@accu-chek.co.za
www.rochediabetescaremea.com

Distributed in the United Kingdom by:

Roche Diabetes Care Limited
Charles Avenue, Burgess Hill
West Sussex, RH15 9RY, **United Kingdom**

Accu-Chek Pump Careline ¹⁾:
UK Freephone number: 0800 731 22 91
ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training
purposes

Some mobile operators may charge for
calls to these numbers.

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